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Power

- Press the power button to switch the device ON and OFF
- When the power is off, press for 3~5 seconds to power on
- When the power is on, press for 3~5 seconds to shut down

Automatic shut-down

- When the device is not connected to the mobile phone, it will automatically shut down after 60 to 70 seconds.

Usage time

- Approximately 48 hours

LED indicator

Device State			LED State	
Non-charging	Power OFF		LED OFF	
	Power ON	BLE not connected	Blinking (3 Hz)	Green light (50% ~ 100%) Yellow light (15% ~ 50%) Red light (0% ~ 15%)
		BLE connected	LED ON	
Charging	Power OFF		Blinking (1 Hz)	
	Power ON	BLE not connected		
		BLE connected		
Updating firmware			LED OFF	

Specification

- Battery: 3.7V, 155mAh lithium battery
- Continuous use: 24 hours
- Charging voltage: 5V DC
- Input impedance: >30M Ohm
- Amplification gain: 303V/V
- Common mode rejection ratio: ≥ 70 dB
- Firmware can be updated by APP

Size

- Device: 70 x 37 x 12mm

Weight

- 27g

Recommended working environment

- Temperature: $-20^{\circ}\text{C} \sim 60^{\circ}\text{C}$
- Humidity: 10% ~ 95%
- Air pressure: 500 ~ 1060 hPa

Precautions

- The device cannot be used with the fast-charging chargers. It is recommended to purchase non-fast-charging USB chargers that comply with IEC 60601-1 or IEC 62368-1 or IEC 60950-1 certification.
- If the battery power is lower than 15%, please charge the device as soon as possible to ensure proper functioning. If the device has not been used for a long time, it needs to be charged at least once every three to six months to maintain the health of the lithium battery.
- After the firmware update is completed, the device will return to the shutdown state. You can restart the device according to the normal startup procedure.

Data transfer mode

- BLE 5.2

ECG and heart rate

- Sampling frequency: 1000Hz
- Resolution: 14-bit
- Measurement range: 15 ~ 300 bpm
- Electrode off detection

Heart Rate Variability (HRV)

- Time domain metrics: SDNN, RMSSD
- Frequency domain metrics: HF, LF, LF/HF, TP